

USP 797 & 800 Compliance Specialists

Everything you need to achieve and maintain a compliant cleanroom environment

Sterile IPA

70% Isopropanol / 30% USP Grade Purified Water



Description

VAI's DECON-AHOL® Products are formulated with USP Isopropyl Alcohol (IPA). DECON-AHOL Products are available in 91%, 70%, and 60% USP IPA with USP Water for Injection (WFI), as well as an aseptically filtered 99% USP IPA. DECON-AHOL WFI® Formula 70% is an EPA registered hard surface disinfectant and sanitizer when used as directed.

DECON-AHOL Products have been designed specifically for pharmaceutical, biotechnology, healthcare, and medical device cleaning rotations for use in both aseptic and non-aseptic environments. DECON-AHOL WFI Formula can be used for the decontamination of numerous surfaces within any cleanroom operation. Each lot of DECON-AHOL is sterility tested according to current USP Compendium, is completely traceable, and has been completely validated for sterility and shelf life. All shipments of DECON-AHOL are delivered with lot specific certificates.

DECON-AHOL is ready to use, available in sterile and non-sterile, and available in a variety of container sizes including aerosol cans, trigger sprays, squeeze bottles, gallon sizes, bulk sizes, and ASEPTI-CLEANSE hands-free dispenser bottles. Our new ASEPTI-CLEANSE Bottles, part numbers DECWFIFOT-01 & DECWFI-BOT-02, CLEANSE Bottles use a new and improved ASEPTI-CLEANSE operating system.

Features and Benefits

- Each sterile container is double bagged
- Quadruple bagged in the ABCD Cleanroom Introduction System®
- Delivered with lot specific Certificate of Analysis, Certificate of Sterility, Certificate of Irradiation, and LAL Test Report
- Individually labeled with lot number and expiration
- Low remaining residue
- Low in endotoxin levels
- · For use on a multitude of surfaces
- Ready-to-use





Product Uses

- Where a sterile alcohol solution that is formulated with USP Water for Injection is required
- In spray and wipe downs on hard, non-porous, inanimate surfaces
- Aseptic filling, gowning rooms, general manufacturing areas, process lines
- Machinery, tools, tables, counters, laminar flow benches, carts, shelves
- Exterior packaging, accessories
- Glass, plastic, vinyl, stainless steel
- Gloves

Quality and Manufacturing

- Raw components are air washed with 0.2 micron filtered air
- Filled in an ISO 5 (Grade A/B, Former Class 100)
- Filtered at 0.2 microns
- Gamma irradiated at a 10-6 SAL
- Lot sterility tested according to current USP compendium
- Completely lot traceable from start to finish
- Completely validated for sterility and shelf life

Product Numbers		
DECWFI-TR-05	32oz Trigger Bottle	12 btls / cs



